


SYSTEM TWO

LED Series

THE PINNACLE OF THE SYSTEM TWO SERIES



- Amazingly cool light
- Ultra high-efficiency – very low electrical consumption
- Excellent color temperature
- 50,000 hour-rated LED life
- 5-stage dimming, on/off switch, focus control and light head positioning – all controlled by our unique sterilizable handle



SYSTEM TWO

LED Series

The **SYSTEM TWO LED** Surgery Lights represent the pinnacle of the System Two series. Utilizing the latest LED technology, these lights are the brightest, longest-lasting light source in our history.

Each light head offers 120,000 lux output with 50,000 hour-rated LED life. The 4,300° K light color offers a beautifully white light field that provides excellent visualization. The five-stage dimming control and 9"– 13" light field are both adjusted entirely with the sterilizable handle, putting all control in the sterile field.

Available in several different mounting configurations, along with video capability, the System Two LED offers a tremendous value in O.R. lighting.

- 120,000 lux output at 1 meter
- 4,300° Kelvin color temperature
- 9"– 13" adjustable pattern
- Five-stage dimming
- 50,000 hour-rated LED life
- Ability to control all functions in the sterile field



INDICATOR PANEL

MONITOR ARM

Our video monitor arm is a convenient solution for monitor placement.

- Sterile handle for adjustments within the sterile field
- Provides options for High Definition, fiber optics cabling, video source switching and compatibility with PACS
- Can be ordered with or without monitor
- Ceiling mount version as solo arm or shared with lights



HIGH-DEFINITION VIDEO CAMERA

Video Camera capability at a very low price.

- 1080P high-definition
- HD-SDI video output
- 20x optical zoom/12x digital zoom
- Auto and manual focus
- Sterilizable handle – standard
- Mounted on its own independent arm assembly
- Can be ordered in combination with one or two lights on a Duo or Trio ceiling mounting



SYSTEM TWO ARM SYSTEM

Available in Solo, Duo and Trio mounts, these arms provide continuous rotation in three key movement points, and provide exceptional reach and articulation. The arms can be configured with light heads and video cameras or monitor arms (monitor and camera arms do not feature continuous rotation).

SYSTEM TWO LED SOLO

The Solo offers unsurpassed longevity and flexibility in a package that is incredibly affordable. The large diameter output provides excellent shadow control – even with just a single light head. This mounting can also support the P1 video camera or M1 monitor arm.



SYSTEM TWO LED DUO

Providing dual lights for superior shadow control and intensity, or the option of one light and one video accessory (monitor or video camera). Without a doubt, an unparalleled value in the surgery market today.



SYSTEM TWO LED TRIO

The Trio represents the ultimate flexibility in the surgery room. Configurations include three light heads, two light heads with video monitor arm, or video camera with two light heads. Complete surgery lighting and video capability with superb performance.



SURGERY LIGHTING

**Unparalleled
Performance**



High Efficiency



Stunning Aesthetics



**Wide Array of
Mounting & Video
Options**



Exceptional Value



Actual photograph of System Two LED light output.

*Medical Illumination International products are sold through medical equipment distributors, worldwide.
Please contact us for referral to a distributor in your region.*

**MEDICAL
ILLUMINATION**
INTERNATIONAL

WHERE ILLUMINATION IS DEFINED

547 Library Street San Fernando, CA 91340 USA
(818) 838-3025 (800) 831-1222 (U.S.A.)
fax (818) 838-3725
www.medillum.com

