



The Colon Hydromat standard shows up in a new design, new color and established quality.

Clear composition of the temperature pre-selection and water flow controller, allows the user a comfortable handling.



Highlights of this system:

- highest safety thanks to proved quality and most modern technology
- precise notification about treatment status on a display
- Temperature pre-selection and water flow controller allows individual settings
- modern Design
- user-friendly control elements guarantee an uncomplicated treatment

Technical data

Name of appliance:	Colon Hydromat standard
Power supply:	100-240 V~50-60 Hz
Water pressure (cold and warm):	2 - 4 bar
Dimensions (BxHxT):	495 x 480 x 250 mm
Weight:	18,2 kg
Product category acc. to MP-directives 93/42 EWG:	II a
	Art.-No. 15570 left Art.-No. 15571 right Art.-No. 15572 left Art.-No. 15573 right Art.-No. 15574 left Art.-No. 15575 right