



Medozon ip

- Developed for intraperitoneal insufflation with ozone
- Simple analysis thanks to internal data storage
- A touch screen and integrated operational guidance guarantee for easy handling
- Multilingual operational guidance (German, English, Spanish)



Technical data

Name:	Medozon ip
Power supply:	230V~/50 Hz or 115V~/50 Hz
Input:	135 VA
Concentration ozone:	5-60 µg/ml (regulation via analyser)
Flow through:	0.8 l/min
Size (WxHxD):	535 x 180 x 460 mm
Weight:	15.5 kg
Classification according to Directive 93/42 EEC	II b

■ Ord. no. 23000
■ Ord. no. 23001

Medozon^{ip}