

SYSTEM TWO

HALOGEN Series



- Excellent Value
- Halogen Light Source
- Four-Stage Dimming
- Multiple Mounting Options
- Video Integration Capability



SYSTEM TWO

HALOGEN Series

The **SYSTEM TWO** Halogen Surgery Lights incorporate proven halogen reliability and an unparalleled list of standard features, making it an industry standard. Full focus control, 4 stage dimming, high color temperature, a complete array of mounting options and low cost make the Halogen Surgery Lights a truly great value.



S2

- 102,000 lux (9500 ft candles) output at 1 meter
- 4,000° Kelvin color temperature

Two light heads are available – the S2, with 20" reflector and the S1 "satellite," with 16" reflector.

Along with the Surgery Light, we offer advanced video integration capability with video camera and video monitoring options.

The **SYSTEM TWO** ceiling mounting and extension arm have been totally redesigned for far easier mounting and improved stability.



S1

- 72,000 lux (6700 ft candles) output at 1 meter
- 4,000° Kelvin color temperature



- 27x optical zoom
- Auto and manual focus
- Optional 8 video input controller (at additional cost)
- Sterilizable handle – standard

P1 VIDEO CAMERA

The **SYSTEM TWO** Video Camera provides excellent video imaging capabilities, at an incredibly low price. Additionally, as the camera is mounted on its own independent arm assembly, it can easily be ordered in combination with one or two Halogen Surgery Lights, creating a Duo or Trio system.

Single camera controller is supplied. An optional 8 input controller is also available at extra cost.



M1 MONITOR ARM

SYSTEM TWO's Video Monitor Arm ideally couples **SYSTEM TWO** Surgery Lights with advanced video imaging, to address the ever-growing needs in the surgery room.

The Monitor Arm can be ordered in combination with ceiling mounted Surgery Lights, as well as a Wall-mounted model, and can be ordered with or without video monitor.



SOLO MOUNT




S2 Pictured

The **SYSTEMTWO** Solo offers great performance and high value, with no true competitor in the market. With the **SYSTEMTWO** articulating arm assembly, which features continuous rotation of all key pivot points, the Solo is the ideal single Surgery Light.

Available in S2 or S1 models.

DUO MOUNT




S2 & S2 Pictured

Our Duo model **SYSTEMTWO** allows for the pairing of two Surgery Lights, the combination of Surgery Light and monitor arm, or Surgery Light and Video Camera. Each articulating arm rotates fully independently of the other, allowing limitless positioning and flexibility.

Both S2 & S1 light heads are available in any of these combinations.



TRIO MOUNT



S2, P1 & S1 Pictured

With our **SYSTEMTWO** Trio, the widest range of capabilities exists. A three Surgery head model is available as well as combinations featuring two Surgery Lights/Video Camera, two Surgery Lights/Monitor and so on.

The Trio provides complete surgery lighting with full video integration to satisfy the most discerning surgery room requirements.



SYSTEM TWO

HALOGEN Series

PORTABLE SURGERY LIGHT



The Portable model offers great flexibility of movement. As well, the head and arm can easily be folded down completely within the footprint of the cantilevered floor base.

WALL MOUNTS



Wall mountings are available for both the Surgery Light and the Video Monitor Arm, providing full flexibility. When a wall mounting is required.

**MEDICAL
ILLUMINATION**
INTERNATIONAL

WHERE ILLUMINATION IS DEFINED

547 Library Street San Fernando, CA 91340 USA
(818) 838-3025 (800) 831-1222 (U.S.A.)
fax (818) 838-3725
www.medillum.com


Manufactured
In The USA